pParse+ User Guide

Friday, January 10, 2014

Contents

pParse+ User Guide	2
Introduction	2
Recommended Minimum System Requirements.....	2
Usage 1.....	2
Usage 2.....	2
How to Set the Parameters?	3
pParse+ 用户手册.....	4
软件简介.....	4
最小系统配置要求.....	4
使用方法 1.....	4
使用方法 2.....	4
参数文件配置说明.....	5

pParse+ User Guide

Introduction

pParse+ is a software tool for extracting MS\MS spectra from *.RAW format. It inspects the isolation windows corresponding to each MS\MS spectra and exports the most possible precursors. Besides pParse+ recalibrates the mass of monoisotopic peak and could detect co-eluted precursors. Currently, users can set co-eluted as '0' to export only one recalibrated precursor monoisotopic peak for each MS\MS scan.

pParse+ do not need install. You can use it after decompressed it into a folder.

Recommended Minimum System Requirements

- CPU: Pentium III or higher
- Main Memory: 1GB
- Operating System: Windows XP/Vista/7/8
- **Xcalibur 2.1 or higher**

Usage 1

a. Set the parameters. A template of parameter configure file named pParse.para can be found in the folder where pParse+ is decompressed.

b. Save the parameter configure file and open command line prompt in current folder, run pParse+ as follow.

```
pParse.exe pParse.para
```

Usage 2

Now pParse has 21 options and each option has a default value. Users need configure a few of them to run pParse. For example, both

```
pParse.exe -D D:\mydata\First.raw
```

and

```
pParse.exe -D D:\mydata\
```

are legal.

The command line

pParse.exe

will print a detailed user guide about usage of more options.

How to Set the Parameters?

To correctly run pParse+, users need to set a few parameters. Some of the parameters are very important. For example, if **co-eluted** is set as 1, pParse will export all co-eluted precursors.

Usually, users only need to set datapath to run pParse.exe.

[Basic Options]

```
datapath=E:\MannData\  
# both folder or file is legal.
```

[Advanced Option]

```
co-elute=0  
# 0, do not export co-eluted precursors.  
# 1, export co-eluted precursors.  
input_format=raw  
# raw/ms1  
isolation_width=2  
# default=2.  
mars_threshold=-0.68  
trainingset=.\TrainingSet.txt  
ipv_file=.\IPV.txt.
```

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